

TOSHIBA AMERICA MEDICAL SYSTEMS, INC.
2441 Michelle Drive, Tustin, CA 92780
Phone: (714) 730-5000

510(k) SUMMARY

MAY 23 2014

1. SUBMITTER'S NAME:
Toshiba America Medical Systems, Inc.

2. ADDRESS:
2441 Michelle Drive
Tustin, CA 92780-2068

3. ESTABLISHMENT REGISTRATION:
2020563

4. CONTACT PERSON:
Paul Biggins
Director, Regulatory Affairs
(714) 669-7808

5. Date Prepared:
March 21, 2014

6. TRADE NAME(S):
Diagnostic Ultrasound System
Aplio Artida (SSH-880CV), V3.2

7. COMMON NAME:
System, Diagnostic Ultrasound

8. DEVICE CLASSIFICATION:
Class II
Ultrasonic Pulsed Doppler Imaging System – Product Code: 90-IYN [per 21 CFR 892.1550]
Ultrasonic Pulsed Echo Imaging System – Product Code: 90-IYO [per 21 CFR 892.1560]
Diagnostic Ultrasonic Transducer – Product Code: 90-ITX [per 21 CFR 892.1570]

9. PREDICATE DEVICE:

Product	Marketed by	510(k) Number	Clearance Date
Aplio Artida (SSH-880CV), V3.0	Toshiba America Medical Systems	K121577	June 22, 2012

10. REASON FOR SUBMISSION:

Modification of a cleared device

11. DEVICE DESCRIPTION:

The APLIO ARTIDA, Model SSH-880CV is a mobile diagnostic ultrasound system. This is a Track 3 device that employs a wide array of probes including convex, pencil, flat linear array and sector array, with a frequency range of approximately 2.0 MHz to 7.5 MHz. This system supports basic measurements including distance, time, angle, and trace, as well as combinations of some basic measurements.

12. INDICATIONS FOR USE:

The DIAGNOSTIC ULTRASOUND SYSTEM APLIO ARTIDA (Model SSH-880CV) is intended to be used for the following types of studies: cardiac, transesophageal, abdominal and peripheral vascular.

13. SUBSTANTIAL EQUIVALENCE:

This device is substantially equivalent to the Aplio Artida (SSH-880CV), V3.0, 510(k) cleared under K121577, marketed by Toshiba America Medical Systems. The **Aplio Artida (SSH-880CV), V3.2**, functions in a manner similar to and is intended for the same use as the predicate device. The subject device includes modifications to the cleared device which improves upon existing features including contour tracing and expansion of the sector display function for two existing transducers and use of a proprietary 4D render. A comparison table is included in this submission detailing the similarities and differences between the predicate device and the subject device.

14. SAFETY:

The device is designed and manufactured under the Quality System Regulations as outlined in 21 CFR § 820 and ISO 13485 Standards. This device is in conformance with the applicable parts of the IEC60601-1, IEC 60601-1-2, IEC 60601-2-37, IEC 62304 and NEMA UD3 Output Display standards.

15. TESTING

Risk Analysis and verification/validation testing conducted through bench testing are included in this submission which demonstrates that the requirements for the improved/added features have been met.

Software Documentation for a Moderate Level of Concern, per the FDA guidance document, "Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices Document" issued on May 11, 2005, is also included as part of this submission.

Additionally, testing of the modified system was conducted in accordance with the applicable standards published by the International Electrotechnical Commission (IEC) for Medical Devices.

16. CONCLUSION

The modifications incorporated into the **Aplio Artida (SSH-880CV), V3.2** do not change the indications for use or the intended use of the device. Based upon bench testing and successful completion of software validation, application of risk management and design controls, it is concluded that this device is safe and effective for its intended use.



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center – WO66-G609
Silver Spring, MD 20993-0002

May 23, 2014

Toshiba Medical Systems Corporation
% Mr. Paul Biggins
Director, Regulatory Affairs
Toshiba America Medical Systems, Inc.
2441 Michelle Drive
TUSTIN CA 92780

Re: K140729

Trade/Device Name: Aplio Artida (SSH-880CV), V3.2
Regulation Number: 21 CFR 892.1550
Regulation Name: Ultrasonic pulsed doppler imaging system
Regulatory Class: II
Product Code: IYN, IYO, ITX
Dated: March 21, 2014
Received: March 24, 2014

Dear Mr. Biggins:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

This determination of substantial equivalence applies to the following transducers intended for use with the Aplio Artida, V3.2, as described in your premarket notification:

Transducer Model Number

PST-25SX	PST-30BT	PST-30SBT
PST-50BT	PST-65AT	PLT-704SBT
PC-20M	PST-25BT	PVT-375BT
PLT-704AT	PET-508MA	PET-510MB
PET-512MC		

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be

found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638 2041 or (301) 796-7100 or at its Internet address


<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRI's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



for

Janine M. Morris
Director
Division of Radiological Health
Office of In Vitro Diagnostics
and Radiological Health
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K140729

Device Name

Aplio Artida™ (SSH-880CV), V3.2

Indications for Use (Describe)

The DIAGNOSTIC ULTRASOUND SYSTEM APLIO ARTIDA (Model SSH-880CV) is intended to be used for the following types of studies: cardiac, transesophageal, abdominal and peripheral vascular.

Type of Use (Select one or both, as applicable)


☒ Prescription Use (Part 21 CFR 801 Subpart D)

☐ Over-The-Counter Use (21 CFR 801 Subpart C)

PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON A SEPARATE PAGE IF NEEDED.

FOR FDA USE ONLY

Concurrence of Center for Devices and Radiological Health (CDRH) (Signature)



This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

Toshiba America Medical Systems, Inc.

510(k) Premarket Notification
Aplio Artida™(v3.2) SSH-880CV Ultrasound System

System: Aplio Artida (SSH-880CV), V3.2
Transducer: _____

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application Specific (Tracks 3)	Mode of Operation											
	B	M	PWD	CWD	Color Doppler	Combined (Specify) *	THI	Dynamic Flow	Power	TDI	CHI 2D	4D (Realtime 3D)
Ophthalmic												
Fetal												
Abdominal	P	P	P	P	P	P	P	P	P	P	P	
Intra-operative (Specify)												
Intra-operative (Neuro)												
Laparoscopic												
Pediatric												
Small Organ (Specify) (1)	P	P	P		P	P	P	P	P			
Neonatal Cephalic												
Adult Cephalic												
Trans-rectal												
Trans-vaginal												
Trans-urethral												
Trans-esoph. (non-Card.)												
Musculo-skeletal (Conventional)	P	P	P		P	P	P	P	P			
Musculo-skeletal (Superficial)	P	P	P		P	P	P	P	P			
Intravascular												
Other (Specify)												
Cardiac Adult	P	P	P	P	P	P	P	P	P	P	P	P
Cardiac Pediatric	P	P	P	P	P	P	P	P	P	P	P	P
Intravascular (Cardiac)												
Trans-esoph. (Cardiac)	P	P	P	P	P	P	P			P		
Intra-cardiac												
Other (Specify)												
Peripheral vessel	P	P	P	P	P	P	P	P	P			
Other (Specify)												

N = new indication; P = previously cleared by FDA; E = added under this appendix
Previous 510k of the transducer: K121577

*Combined mode includes B/M; B/PWD; BDF/PWD; BDF/MDF; BDF/MDF/PWD; 2D/CWD; BDF/CWD

Toshiba America Medical Systems, Inc.

510(k) Premarket Notification
Aplio Artida™(v3.2) SSH-880CV Ultrasound System

System: Aplio Artida (SSH-880CV), V3.2
Transducer: PST-25SX

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation											
Specific (Tracks 3)	B	M	PWD	CWD	Color Doppler	Combined (Specify) *	THI	Dynamic Flow	Power	TDI	CHI 2D	4D (Realtime 3D)
Ophthalmic												
Fetal												
Abdominal												
Intra-operative (Specify)												
Intra-operative (Neuro)												
Laparoscopic												
Pediatric												
Small Organ (Specify) (1)												
Neonatal Cephalic												
Adult Cephalic												
Trans-rectal												
Trans-vaginal												
Trans-urethral												
Trans-esoph. (non-Card.)												
Musculo-skeletal (Conventional)												
Musculo-skeletal (Superficial)												
Intravascular												
Other (Specify)												
Cardiac Adult	P						P					P
Cardiac Pediatric	P						P					P
Intravascular (Cardiac)												
Trans-esoph. (Cardiac)												
Intra-cardiac												
Other (Specify)												
Peripheral vessel												
Other (Specify)												

N = new indication; P = previously cleared by FDA; E = added under this appendix
Previous 510(k) of the transducer: K121577

*Combined mode includes B/M; B/PWD; BDF/PWD; BDF/MDF; BDF/MDF/PWD; 2D/CWD; BDF/CWD

System: Aplio Artida (SSH-880CV), V3.2V
 Transducer: PST-30BT

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application Specific (Tracks 3)	Mode of Operation											
	B	M	PWD	CWD	Color Doppler	Combined (Specify)	THI	Dynamic Flow	Power	TDI	CHI 2D	4D (Realtime 3D)
Ophthalmic												
Fetal												
Abdominal	P	P	P	P	P	P	P	P	P	P	P	
Intra-operative (Specify)												
Intra-operative (Neuro)												
Laparoscopic												
Pediatric												
Small Organ (Specify) (1)												
Neonatal Cephalic												
Adult Cephalic												
Trans-rectal												
Trans-vaginal												
Trans-urethral												
Trans-esoph. (non-Card.)												
Musculo-skeletal (Conventional)												
Musculo-skeletal (Superficial)												
Intravascular												
Other (Specify)												
Cardiac Adult	P	P	P	P	P	P	P	P	P	P	P	
Cardiac Pediatric	P	P	P	P	P	P	P	P	P	P	P	
Intravascular (Cardiac)												
Trans-esoph. (Cardiac)												
Intra-cardiac												
Other (Specify)												
Peripheral vessel												
Other (Specify)												

N = new indication; P = previously cleared by FDA; E = added under this appendix
 Previous 510(k) of the transducer: K121577

*Combined mode includes B/M; B/PWD; BDF/PWD; BDF/MDF; BDF/MDF/PWD; 2D/CWD; BDF/CWD

System: Aplio Artida (SSH-880CV), V3.2
Transducer: PST-30SBT

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application Specific (Tracks 3)	Mode of Operation											
	B	M	PWD	CWD	Color Doppler	Combined (Specify)	THI	Dynamic Flow	Power	TDI	CHI 2D	4D (Realtime 3D)
Ophthalmic												
Fetal												
Abdominal												
Intra-operative (Specify)												
Intra-operative (Neuro)												
Laparoscopic												
Pediatric												
Small Organ (Specify) (1)												
Neonatal Cephalic												
Adult Cephalic												
Trans-rectal												
Trans-vaginal												
Trans-urethral												
Trans-esoph. (non-Card.)												
Musculo-skeletal (Conventional)												
Musculo-skeletal (Superficial)												
Intravascular												
Other (Specify)												
Cardiac Adult	P	P	P	P	P	P	P			P	P	
Cardiac Pediatric	P	P	P	P	P	P	P			P	P	
Intravascular (Cardiac)												
Trans-esoph. (Cardiac)												
Intra-cardiac												
Other (Specify)												
Peripheral vessel												
Other (Specify)												

N = new indication; P = previously cleared by FDA; E = added under this appendix
Previous 510(k) of the transducer: K121577

*Combined mode includes B/M; B/PWD; BDF/PWD; BDF/MDF; BDF/MDF/PWD; 2D/CWD; BDF/CWD

Toshiba America Medical Systems, Inc.

510(k) Premarket Notification
Aplio Artida™(v3.2) SSH-880CV Ultrasound System

System: Aplio Artida (SSH-880CV), V3.2
Transducer: PST-50BT

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application Specific (Tracks 3)	Mode of Operation											
	B	M	PWD	CWD	Color Doppler	Combined (Specify) *	THI	Dynamic Flow	Power	TDI	CHI 2D	4D (Realtime 3D)
Ophthalmic												
Fetal												
Abdominal												
Intra-operative (Specify)												
Intra-operative (Neuro)												
Laparoscopic												
Pediatric												
Small Organ (Specify) (1)												
Neonatal Cephalic												
Adult Cephalic												
Trans-rectal												
Trans-vaginal												
Trans-urethral												
Trans-esoph. (non-Card.)												
Musculo-skeletal (Conventional)												
Musculo-skeletal (Superficial)												
Intravascular												
Other (Specify)												
Cardiac Adult	P	P	P	P	P	P	P			P		
Cardiac Pediatric	P	P	P	P	P	P	P			P		
Intravascular (Cardiac)												
Trans-esoph. (Cardiac)												
Intra-cardiac												
Other (Specify)												
Peripheral vessel												
Other (Specify)												

N = new indication; P = previously cleared by FDA; E = added under this appendix
Previous 510(k) of the transducer: K121577

*Combined mode includes B/M; B/PWD; BDF/PWD; BDF/MDF; BDF/MDF/PWD; 2D/CWD; BDF/CWD

Toshiba America Medical Systems, Inc.

510(k) Premarket Notification
Aplio Artida™(v3.2) SSH-880CV Ultrasound System

System: Aplio Artida (SSH-880CV), V3.2
Transducer: PST-65AT

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation											
Specific (Tracks 3)	B	M	PWD	CWD	Color Doppler	Combined (Specify) *	THI	Dynamic Flow	Power	TDI	CHI 2D	4D (Realtime 3D)
Ophthalmic												
Fetal												
Abdominal												
Intra-operative (Specify)												
Intra-operative (Neuro)												
Laparoscopic												
Pediatric												
Small Organ (Specify) (I)												
Neonatal Cephalic												
Adult Cephalic												
Trans-rectal												
Trans-vaginal												
Trans-urethral												
Trans-esoph. (non-Card.)												
Musculo-skeletal (Conventional)												
Musculo-skeletal (Superficial)												
Intravascular												
Other (Specify)												
Cardiac Adult	P	P	P	P	P	P	P			P		
Cardiac Pediatric	P	P	P	P	P	P	P			P		
Intravascular (Cardiac)												
Trans-esoph. (Cardiac)												
Intra-cardiac												
Other (Specify)												
Peripheral vessel												
Other (Specify)												

N = new indication; P = previously cleared by FDA; E = added under this appendix
Previous 510(k) of the transducer: K121577

*Combined mode includes B/M; B/PWD; BDF/PWD; BDF/MDF; BDF/MDF/PWD; 2D/CWD; BDF/CWD

Toshiba America Medical Systems, Inc.

510(k) Premarket Notification
Aplio Artida™(v3.2) SSH-880CV Ultrasound System

System: Aplio Artida (SSH-880CV), V3.2
Transducer: PLT-704SBT

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application Specific (Tracks 3)	Mode of Operation											
	B	M	PWD	CWD	Color Doppler	Combined (Specify)	THI	Dynamic Flow	Power	TDI	CHI 2D	4D (Realtime 3D)
Ophthalmic												
Fetal												
Abdominal												
Intra-operative (Specify)												
Intra-operative (Neuro)												
Laparoscopic												
Pediatric												
Small Organ (Specify) (1)	P	P	P		P	P	P	P	P			
Neonatal Cephalic												
Adult Cephalic												
Trans-rectal												
Trans-vaginal												
Trans-urethral												
Trans-esoph. (non-Card.)												
Musculo-skeletal (Conventional)	P	P	P		P	P	P	P	P			
Musculo-skeletal (Superficial)	P	P	P		P	P	P	P	P			
Intravascular												
Other (Specify)												
Cardiac Adult												
Cardiac Pediatric												
Intravascular (Cardiac)												
Trans-esoph. (Cardiac)												
Intra-cardiac												
Other (Specify)												
Peripheral vessel	P	P	P		P	P	P	P	P			
Other (Specify)												

N = new indication; P = previously cleared by FDA; E = added under this appendix
Previous 510(k) of the transducer: K121577

*Combined mode includes B/M; B/PWD; BDF/PWD; BDF/MDF; BDF/MDF/PWD; 2D/CWD; BDF/CWD

System: Aplio Artida (SSH-880CV), V3.2
 Transducer: PC-20M

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application Specific (Tracks 3)	Mode of Operation											
	B	M	PWD	CWD	Color Doppler	Combined (Specify) *	THI	Dynamic Flow	Power	TDI	CHI 2D	4D (Realtime 3D)
Ophthalmic												
Fetal												
Abdominal												
Intra-operative (Specify)												
Intra-operative (Neuro)												
Laparoscopic												
Pediatric												
Small Organ (Specify) (1)												
Neonatal Cephalic												
Adult Cephalic												
Trans-rectal												
Trans-vaginal												
Trans-urethral												
Trans-esoph. (non-Card.)												
Musculo-skeletal (Conventional)												
Musculo-skeletal (Superficial)												
Intravascular												
Other (Specify)												
Cardiac Adult				P								
Cardiac Pediatric				P								
Intravascular (Cardiac)												
Trans-esoph. (Cardiac)												
Intra-cardiac												
Other (Specify)												
Peripheral vessel				P								
Other (Specify)												

N = new indication; P = previously cleared by FDA; E = added under this appendix
 Previous 510(k) of the transducer: K121577

*Combined mode includes B/M; B/PWD; BDF/PWD; BDF/MDF; BDF/MDF/PWD; 2D/CWD; BDF/CWD

System: Aplio Artida (SSH-880CV), V3.2
 Transducer: PST-25BT

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation											
Specific (Tracks 3)	B	M	PWD	CWD	Color Doppler	Combined (Specify) *	THI	Dynamic Flow	Power	TDI	CHI 2D	4D (Realtime 3D)
Ophthalmic												
Fetal												
Abdominal												
Intra-operative (Specify)												
Intra-operative (Neuro)												
Laparoscopic												
Pediatric												
Small Organ (Specify) (1)												
Neonatal Cephalic												
Adult Cephalic												
Trans-rectal												
Trans-vaginal												
Trans-urethral												
Trans-esoph. (non-Card.)												
Musculo-skeletal (Conventional)												
Musculo-skeletal (Superficial)												
Intravascular												
Other (Specify)												
Cardiac Adult	P	P	P	P	P	P	P			P	P	
Cardiac Pediatric	P	P	P	P	P	P	P			P	P	
Intravascular (Cardiac)												
Trans-esoph. (Cardiac)												
Intra-cardiac												
Other (Specify)												
Peripheral vessel												
Other (Specify)												

N = new indication; P = previously cleared by FDA; E = added under this appendix
 Previous 510(k) of the transducer: K121577

*Combined mode includes B/M; B/PWD; BDF/PWD; BDF/MDF; BDF/MDF/PWD; 2D/CWD; BDF/CWD

System: Aplio Artida (SSH-880CV), V3.2
 Transducer: PVT-375BT

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application Specific (Tracks 3)	Mode of Operation											
	B	M	PWD	CWD	Color Doppler	Combined (Specify) *	THI	Dynamic Flow	Power	TDI	CHI 2D	4D (Realtime 3D)
Ophthalmic												
Fetal												
Abdominal	P	P	P		P	P	P	P	P			
Intra-operative (Specify)												
Intra-operative (Neuro)												
Laparoscopic												
Pediatric												
Small Organ (Specify) (1)												
Neonatal Cephalic												
Adult Cephalic												
Trans-rectal												
Trans-vaginal												
Trans-urethral												
Trans-esoph. (non-Card.)												
Musculo-skeletal (Conventional)												
Musculo-skeletal (Superficial)												
Intravascular												
Other (Specify)												
Cardiac Adult												
Cardiac Pediatric												
Intravascular (Cardiac)												
Trans-esoph. (Cardiac)												
Intra-cardiac												
Other (Specify)												
Peripheral vessel	P	P	P		P	P	P	P	P			
Other (Specify)												

N = new indication; P = previously cleared by FDA; E = added under this appendix
 Previous 510(k) of the transducer: K121577

*Combined mode includes B/M; B/PWD; BDF/PWD; BDF/MDF; BDF/MDF/PWD; 2D/CWD; BDF/CWD

Toshiba America Medical Systems, Inc.

510(k) Premarket Notification
Aplio Artida™(v3.2) SSH-880CV Ultrasound System

System: Aplio Artida (SSH-880CV), V3.2
Transducer: PLT-704AT

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation											
Specific (Tracks 3)	B	M	PWD	CWD	Color Doppler	Combined (Specify)	THI	Dynamic Flow	Power	TDI	CHI 2D	4D (Realtime 3D)
Ophthalmic												
Fetal												
Abdominal												
Intra-operative (Specify)												
Intra-operative (Neuro)												
Laparoscopic												
Pediatric												
Small Organ (Specify) (1)												
Neonatal Cephalic												
Adult Cephalic												
Trans-rectal												
Trans-vaginal												
Trans-urethral												
Trans-esoph. (non-Card.)												
Musculo-skeletal (Conventional)												
Musculo-skeletal (Superficial)												
Intravascular												
Other (Specify)												
Cardiac Adult												
Cardiac Pediatric												
Intravascular (Cardiac)												
Trans-esoph. (Cardiac)												
Intra-cardiac												
Other (Specify)												
Peripheral vessel	P	P	P		P	P	P	P	P			
Other (Specify)												

N = new indication; P = previously cleared by FDA; E = added under this appendix
Previous 510(k) of the transducer: K121577

*Combined mode includes B/M; B/PWD; BDF/PWD; BDF/MDF; BDF/MDF/PWD; 2D/CWD; BDF/CWD

Toshiba America Medical Systems, Inc.

510(k) Premarket Notification
Aplio Artida™(v3.2) SSH-880CV Ultrasound System

System: Aplio Artida (SSH-880CV), V3.2
Transducer: PET-508MA

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation											
Specific (Tracks 3)	B	M	PWD	CWD	Color Doppler	Combined (Specify)	THI	Dynamic Flow	Power	TDI	CHI 2D	4D (Realtime 3D)
Ophthalmic												
Fetal												
Abdominal												
Intra-operative (Specify)												
Intra-operative (Neuro)												
Laparoscopic												
Pediatric												
Small Organ (Specify) (1)												
Neonatal Cephalic												
Adult Cephalic												
Trans-rectal												
Trans-vaginal												
Trans-urethral												
Trans-esoph. (non-Card.)												
Musculo-skeletal (Conventional)												
Musculo-skeletal (Superficial)												
Intravascular												
Other (Specify)												
Cardiac Adult												
Cardiac Pediatric												
Intravascular (Cardiac)												
Trans-esoph. (Cardiac)	P	P	P	P	P	P	P			P		
Intra-cardiac												
Other (Specify)												
Peripheral vessel												
Other (Specify)												

N = new indication; P = previously cleared by FDA; E = added under this appendix
Previous 510(k) of the transducer: K121577

*Combined mode includes B/M; B/PWD; BDF/PWD; BDF/MDF; BDF/MDF/PWD; 2D/CWD; BDF/CWD

Toshiba America Medical Systems, Inc.

510(k) Premarket Notification
 Aplio Artida™(v3.2) SSH-880CV Ultrasound System

System: Aplio Artida (SSH-880CV), V3.2
 Transducer: PET-S10MB

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation											
Specific (Tracks 3)	B	M	PWD	CWD	Color Doppler	Combined (Specify) *	THI	Dynamic Flow	Power	TDI	CHI 2D	4D (Realtime 3D)
Ophthalmic												
Fetal												
Abdominal												
Intra-operative (Specify)												
Intra-operative (Neuro)												
Laparoscopic												
Pediatric												
Small Organ (Specify) (1)												
Neonatal Cephalic												
Adult Cephalic												
Trans-rectal												
Trans-vaginal												
Trans-urethral												
Trans-esoph. (non-Card.)												
Musculo-skeletal (Conventional)												
Musculo-skeletal (Superficial)												
Intravascular												
Other (Specify)												
Cardiac Adult												
Cardiac Pediatric												
Intravascular (Cardiac)												
Trans-esoph. (Cardiac)	P	P	P	P	P	P	P			P		
Intra-cardiac												
Other (Specify)												
Peripheral vessel												
Other (Specify)												

N = new indication; P = previously cleared by FDA; E = added under this appendix
 Previous 510(k) of the transducer: K121577

*Combined mode includes B/M; B/PWD; BDF/PWD; BDF/MDF; BDF/MDF/PWD; 2D/CWD; BDF/CWD

Toshiba America Medical Systems, Inc.

510(k) Premarket Notification
Aplio Artida™(v3.2) SSH-880CV Ultrasound System

System: Aplio Artida (SSH-880CV), V3.2
Transducer: PET-512MC

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application Specific (Tracks 3)	Mode of Operation											
	B	M	PWD	CWD	Color Doppler	Combined (Specify) *	THI	Dynamic Flow	Power	TDI	CHI 2D	4D (Realtime 3D)
Ophthalmic												
Fetal												
Abdominal												
Intra-operative (Specify)												
Intra-operative (Neuro)												
Laparoscopic												
Pediatric												
Small Organ (Specify) (1)												
Neonatal Cephalic												
Adult Cephalic												
Trans-rectal												
Trans-vaginal												
Trans-urethral												
Trans-esoph. (non-Card.)												
Musculo-skeletal (Conventional)												
Musculo-skeletal (Superficial)												
Intravascular												
Other (Specify)												
Cardiac Adult												
Cardiac Pediatric												
Intravascular (Cardiac)												
Trans-esoph. (Cardiac)	P	P	P	P	P	P	P			P		
Intra-cardiac												
Other (Specify)												
Peripheral vessel												
Other (Specify)												

N = new indication; P = previously cleared by FDA; E = added under this appendix
Previous 510(k) of the transducer: K121577

*Combined mode includes B/M; B/PWD; BDF/PWD; BDF/MDF; BDF/MDF/PWD; 2D/CWD; BDF/CWD